

Infrared micro sterilizer



Infrared micro sterilizer

Product codes:

Reference: B1000

Product description:

Infrared micro sterilizer

Infrared heating to produce 815 °C or 1500° F

Complete infrared sterilization in 5/7 seconds

Safe, convenient and no spatter.

Ideal for platinum loops, needles and tube mouths

Product features:

Dimensions: Well 20 cm deep x 5 cm diameter

Heating up time: 20 min

Working temperature: 260° C (approximately 500° F)

Time: Sterilization : 2-5 seconds

Available models

